




up date 03.03.2015

Super Low Cost Simple DNA extraction protocol for PCR (ver1.9)

1. Place approx. a punch leaf (1/4" hole diameter) into 1.5ml tube containing 100 μ L KAZU buffer (4 drops)



*This can be stored at room temperature for months.

2. Crush leaf disc with blue pestle (best) or macerate with pipette tip (easy).



3. Transfer 5 μ L or 10 μ L supernatant to new tube (avoid chunks of tissue).
4. Add 1ml 95% ethanol, invert tube a few times.
5. Spin down DNA pellet at 12,000 rpm for 30 sec.
6. Discard supernatant.
7. Dry pellet on bench.
8. Dissolve DNA pellet in H₂O (50 μ L or 100 μ L).

PCR

DNA 2.5ul

A primer 2.5ul

B primer 2.5ul

Go taq mix (2X) 7.5ul

95°C 3min-(95°C 30sec-55°C 30sec-72°C 1min) x35-72°C 5min -12°C

Any feedback would be appreciated!

Kazuhiro Kikuchi
Research Scientist
Thomas Brutnell Laboratory
Donald Danforth Plant Science Center
975 N. Warson Rd
St. Louis, MO 63132
Phone: (314) 587-1276
Fax: (314) 587-1376
Email: kkikuchi@danforthcenter.org